

# TMED4 antibody



Catalog Number: 116087

## Product name

TMED4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TMED4 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 26 kDa

## Preparation

This antibody was obtained by immunization of TMED4 recombinant protein (Accession Number: NM\_182547). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116087(TMED4 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:116087(TMED4 Antibody) at dilution of 1:25 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)