

# TMEM108 antibody



Catalog Number: 116092

## Product name

TMEM108 antibody

## Specificity

Human; other species not tested.

## Antibody description

TMEM108 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 70-75 kDa

## Preparation

This antibody was obtained by immunization of TMEM108 recombinant protein (Accession Number: XM\_011513097). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116092(TMEM108 antibody) at dilution of 1:200