

Thyroglobulin antibody

Catalog Number: 116117

Product name

Thyroglobulin antibody

Specificity

Human; other species not tested.

Antibody description

Thyroglobulin Rabbit Polyclonal antibody. Positive IHC detected in human thyroid tissue, human thyroid cancer tissue. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of Thyroglobulin recombinant protein (Accession Number: NM_003235). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


Validations



Immunohistochemistry of paraffin-embedded human thyroid tissue slide using Catalog No:116117(TG Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human thyroid tissue slide using Catalog No:116117(TG Antibody) at dilution of 1:200 (under 40x lens).



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116117(TG Antibody) at dilution of 1:1000