

TIMM17A antibody



Catalog Number: 116131

Product name

TIMM17A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TIMM17A Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer.

Preparation

This antibody was obtained by immunization of TIMM17A recombinant protein (Accession Number: NM_006335). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116131(TIMM17A antibody) at dilution of 1:100 (under 10x lens)