

TIMMDC1 antibody



Catalog Number: 116137

Product name

TIMMDC1 antibody

Specificity

Human; other species not tested.

Antibody description

TIMMDC1 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Observed molecular weight by Western-blot: 32-35 kDa

Preparation

This antibody was obtained by immunization of TIMMDC1 recombinant protein (Accession Number: NM_016589). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:116137(TIMMDC1 antibody) at dilution of 1:1000