

TMEM208 antibody



Catalog Number: 116176

Product name

TMEM208 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TMEM208 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 19-29 kDa

Preparation

This antibody was obtained by immunization of TMEM208 recombinant protein (Accession Number: NM_014187). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116176(TMEM208 Antibody) at dilution of 1:300