

TMEM214 antibody

Catalog Number: 116178

Product name

TMEM214 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMEM214 Rabbit Polyclonal antibody. Positive IF detected in human muscle slides. Positive WB detected in C2C12. Observed molecular weight by Western-blot: 77 kDa

Preparation

This antibody was obtained by immunization of TMEM214 recombinant protein (Accession Number: BC002467). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF, WB

Dilutions


Recommended Dilution:

WB: N/AIF: N/A

Validations



WB result of TMEM214 antibody (Catalog No:116178) with C2C12 cells (Control and KD).



Human Muscle Section

IF result of TMEM214 antibody (Catalog No:116178) with human muscle slides.