

# TMEM50B antibody



Catalog Number: 116190

---

## Product name

TMEM50B antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TMEM50B Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 14kd

## Preparation

This antibody was obtained by immunization of TMEM50B recombinant protein (Accession Number: XM\_011529746). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116190(TMEM50B antibody) at dilution of 1:200