

# TMEM70 antibody



Catalog Number: 116195

## Product name

TMEM70 antibody

## Specificity

Human; other species not tested.

## Antibody description

TMEM70 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive WB detected in A2780 cells. Observed molecular weight by Western-blot: 18 kDa

## Preparation

This antibody was obtained by immunization of TMEM70 recombinant protein (Accession Number: NM\_017866). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:116195(TMEM70 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human heart using Catalog No:116195(TMEM70 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:116195(TMEM70 antibody) at dilution of 1:100 (under 40x lens)