

TPPP3 antibody

Catalog Number: 116214

Product name

TPPP3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TPPP3 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 19-20 kDa

Preparation

This antibody was obtained by immunization of TPPP3 recombinant protein (Accession Number: NM_016140). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



WB result of anti-TPPP3 (Catalog No:116214, 1:500) with HEK-293 cells (siRNA transfected with) by Dr. Yupeng Chen. (Negative control or TPPP3 specific siRNAs were transfected into HEK-293 cells with the final concentration of 20 nM. Forty eight hours after transfection, cells were lysed with RIPA buffer and subjected to SDS



PAGE)

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116214(TPPP3 antibody) at dilution of 1:100