

# TPRKB antibody



Catalog Number: 116215

---

## Product name

TPRKB antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TPRKB Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 20kd

## Preparation

This antibody was obtained by immunization of TPRKB recombinant protein (Accession Number: NM\_001330389). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116215(TPRKB antibody) at dilution of 1:300