

TPST2 antibody



Catalog Number: 116219

Product name

TPST2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TPST2 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of TPST2 recombinant protein (Accession Number: NM_003595). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116219(TPST2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116219(TPST2 antibody) at dilution of 1:100 (under 10x lens)