

TRAF5 antibody



Catalog Number: 116233

Product name

TRAF5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TRAF5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse pancreas tissue. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of TRAF5 recombinant protein (Accession Number: NM_145759). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116233(TRAF5 Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116233(TRAF5 antibody) at dilution of 1:100 (under 10x lens)