

TOP1MT antibody



Catalog Number: 116245

Product name

TOP1MT antibody

Specificity

Human; other species not tested.

Antibody description

TOP1MT Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of TOP1MT recombinant protein (Accession Number: NM_052963). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116245(TOP1MT Antibody) at dilution of 1:1000