

# TOR3A antibody



Catalog Number: 116254

## Product name

TOR3A antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

TOR3A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 46kd

## Preparation

This antibody was obtained by immunization of TOR3A recombinant protein (Accession Number: NM\_022371). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116254(TOR3A antibody) at dilution of 1:300