

TPH2 antibody

Catalog Number: 116264

Product name

TPH2 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

TPH2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of TPH2 recombinant protein (Accession Number: BC114442). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116264(TPH2 Antibody) at dilution of 1:300