

TMPRSS5 antibody



Catalog Number: 116277

Product name

TMPRSS5 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

TMPRSS5 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of TMPRSS5 recombinant protein (Accession Number: BC121803). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116277(TMPRSS5 Antibody) at dilution of 1:1000