

TMPRSS6 antibody



Catalog Number: 116278

Product name

TMPRSS6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMPRSS6 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells, MCF7 cells, mouse liver tissue. Observed molecular weight by Western-blot: 52kd

Preparation

This antibody was obtained by immunization of TMPRSS6 recombinant protein (Accession Number: NM_153019). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116278(TMPRSS6 antibody) at dilution of 1:500