

TMSB4X antibody



Catalog Number: 116280

Product name

TMSB4X antibody

Specificity

Human; other species not tested.

Antibody description

TMSB4X Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 5 kDa

Preparation

This antibody was obtained by immunization of TMSB4X recombinant protein (Accession Number: BC001631). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116280(TMSB4X antibody) at dilution of 1:2000