

TMX4 antibody



Catalog Number: 116284

Product name

TMX4 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

TMX4 Rabbit Polyclonal antibody. Positive FC detected in HepG2 cells. Positive IHC detected in human brain tissue, human testis tissue. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue, A375 cells, human heart tissue, L02 cells. Observed molecular weight by Western-blot: 40-45kd

Validations




Preparation

This antibody was obtained by immunization of TMX4 recombinant protein (Accession Number: NM_021156). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116284(TMx4 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:116284(TMx4 antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF, FC

Dilutions


Recommended Dilution:



TMX4 antibody


Catalog Number: 116284

Immunohistochemical of paraffin-embedded human brain using Catalog No:116284(TM_X4 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:116284(TM_X4 Antibody) at dilution of 1:50 and and Rhodamine-labeled goat anti-rabbit

IgG (red).



1X10⁶ HepG2 cells were stained with .2ug TMX4 antibody (Catalog No:116284, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.