

# TNF-alpha antibody

Catalog Number: 116288

## Product name

TNF-alpha antibody

## Specificity

Human; other species not tested.

## Antibody description

TNF-alpha Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human liver cancer tissue, human tonsillitis tissue.

## Preparation

This antibody was obtained by immunization of TNF-alpha recombinant protein (Accession Number: NM\_000594). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116288(TNF-a Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116288(TNF-a Antibody) at dilution of 1:50 (under 40x lens)