

TNFRSF6B antibody

Catalog Number: 116292

Product name

TNFRSF6B antibody

Specificity

Human; other species not tested.

Antibody description

TNFRSF6B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of TNFRSF6B recombinant protein (Accession Number: NM_003823). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-TNFRSF6B (IP:Catalog No:116292, 4ug; Detection:Catalog No:116292 1:300) with HeLa cells lysate 2480ug.