

# TRIM15 antibody



Catalog Number: 116299

## Product name

TRIM15 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TRIM15 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 52kd

## Preparation

This antibody was obtained by immunization of TRIM15 recombinant protein (Accession Number: NM\_033229). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116299(TRIM15 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using TRIM15 antibody Catalog No:116299 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).