

# TRIM31 antibody



Catalog Number: 116310

## Product name

TRIM31 antibody

## Immunogen

[Human TRIM31 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

TRIM31 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive IP detected in COLO 320 cells. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of TRIM31 recombinant protein (Accession Number: NM\_007028). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:1000

IHC: 1:20-1:200

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116310(TRIM31 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116310(TRIM31 antibody) at dilution of 1:50

# TRIM31 antibody

Catalog Number: 116310

---

IP Result of anti-TRIM31 (IP:Catalog No:116310, 4ug; Detection:Catalog No:116310 1:400) with COLO 320 cells lysate 3200ug.

