

# TRIP13 antibody



Catalog Number: 116327

## Product name

TRIP13 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TRIP13 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HT-1080 cells, human testis tissue, mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of TRIP13 recombinant protein (Accession Number: XM\_011514163). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116327 (TRIP13 antibody) at dilution of 1:400



IP Result of anti-TRIP13 (IP:Catalog No:116327, 3ug; Detection:Catalog No:116327 1:400) with mouse testis tissue lysate 6000ug.