

# TRAPPC9 antibody



Catalog Number: 116346

## Product name

TRAPPC9 antibody

Recommended Dilution:

WB: 1:200-1:2000

## Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

## Antibody description

TRAPPC9 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human heart tissue. Positive IF detected in HEK-293 cells. Positive WB detected in mouse brain tissue, human brain tissue, human skeletal muscle tissue, mouse skeletal muscle tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 128-140 kDa

## Validations




## Preparation

This antibody was obtained by immunization of TRAPPC9 recombinant protein (Accession Number: BC006206). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116346(TRAPPC9,NIBP antibody) at dilution of 1:800

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human kidney using Catalog No:116346(TRAPPC9,NIBP antibody) at dilution of 1:50 (under 10x lens)


## Applications

ELISA, WB, IHC, IF, IP

## Dilutions


# TRAPPC9 antibody

Catalog Number: 116346




Immunohistochemical of paraffin-embedded human kidney using Catalog No:116346(TRAPPC9,NIBP antibody) at dilution of

IP Result of anti-TRAPPC9,NIBP (IP:Catalog No:116346, 4ug; Detection:Catalog No:116346 1:500) with mouse brain tissue lysate 3600ug.



Immunofluorescent analysis of HEK-293 cells using Catalog No:116346(TRAPPC9,NIBP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



1:50 (under 40x lens)