

# TRBV5-4 antibody



Catalog Number: 116353

## Product name

TRBV5-4 antibody

## Specificity

human; other species not tested.

## Antibody description

TRBV5-4 Rabbit Polyclonal antibody. Positive WB detected in U-937 cells. Observed molecular weight by Western-blot: 30-34kd

## Preparation

This antibody was obtained by immunization of TRBV5-4 recombinant protein (Accession Number: BC028083). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:116353(TRBV5-4 antibody) at dilution of 1:400