

# TREX1 antibody

Catalog Number: 116359

## Product name

TREX1 antibody

## Specificity

Human; other species not tested.

## Antibody description

TREX1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of TREX1 recombinant protein (Accession Number: NM\_016381). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116359(TREX1



Antibody) at dilution of 1:300

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116359(TREX1 Antibody) at dilution of 1:300



IP Result of anti-TREX1 (IP:Catalog No:116359, 4ug; Detection:Catalog No:116359 1:300) with HeLa cells lysate 1200ug.