

TSPAN33 antibody

Catalog Number: 116362

Product name

TSPAN33 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSPAN33 Rabbit Polyclonal antibody. Positive IHC detected in human hepatocirrhosis tissue. Positive WB detected in human brain tissue, human blood tissue. Observed molecular weight by Western-blot: 64-72kd

Preparation

This antibody was obtained by immunization of TSPAN33 recombinant protein (Accession Number: NM_178562). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116362(TSPAN33 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116362(TSPAN33 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog

TSPAN33 antibody



Catalog Number: 116362

No:116362(TSPAN33 antibody) at dilution of 1:50 (under 40x lens)