

# TRIP6 antibody



Catalog Number: 116391

## Product name

TRIP6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TRIP6 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HeLa cells, HepG2 cells, human testis tissue, mouse lung tissue. Observed molecular weight by Western-blot: 50-55 kDa

## Preparation

This antibody was obtained by immunization of TRIP6 recombinant protein (Accession Number: NM\_003302). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116391 (TRIP6 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using TRIP6 antibody Catalog No:116391 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).