

TrkC antibody

Catalog Number: 116394

Product name

TrkC antibody

Specificity

Human; other species not tested.

Antibody description

TrkC Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 68kd and 94kd and 145kd

Preparation

This antibody was obtained by immunization of TrkC recombinant protein (Accession Number: XM_011521637). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116394(NTRK3 antibody) at dilution of 1:500