

TSGA13 antibody



Catalog Number: 116431

Product name

TSGA13 antibody

Specificity

Human; other species not tested.

Antibody description

TSGA13 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of TSGA13 recombinant protein (Accession Number: NM_052933). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:116431(TSGA13 Antibody) at dilution of 1:50 (under 10x lens)