

TSN antibody

Catalog Number: 116436

Product name

TSN antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TSN Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in Jurkat cells, HeLa cells, MCF-7 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004622). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:1000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116436(TSN Antibody) at dilution of 1:600



IP Result of anti-TSN (IP:Catalog No:116436, 3ug; Detection:Catalog No:116436 1:300) with MCF-7 cells lysate 1600ug.



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:116436(TSN Antibody) at dilution of 1:50 (under 40x lens)