

TSP50 antibody



Catalog Number: 116439

Product name

TSP50 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse; other species not tested.

IF: 1:10-1:100

Antibody description

TSP50 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HEK-293 cells, human testis tissue, mouse testis tissue. Positive IHC detected in human ovary tumor tissue, human breast cancer tissue, human testis tissue. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 50-55 kDa

Validations




Preparation

This antibody was obtained by immunization of TSP50 recombinant protein (Accession Number: NM_013270). Purification method: Antigen affinity purified.

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116439(TSP50 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116439(TSP50 antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:



TSP50 antibody



Catalog Number: 116439

Immunofluorescent analysis of MCF-7 cells using 1:25 and Alexa Fluor 488-conjugated AffiniPure
Catalog No:116439(TSP50 Antibody) at dilution of Goat Anti-Rabbit IgG(H+L)