

# TSPAN10 antibody

Catalog Number: 116441

## Product name

TSPAN10 antibody

## Immunogen

[Human TSPAN10 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

TSPAN10 Rabbit Polyclonal antibody. Positive IF detected in hTERT-RPE1 cells. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of TSPAN10 recombinant protein (Accession Number: NM\_031945). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



IF result (nuclear and the base of cilia stain) of anti-TSPAN10 (Catalog No:116441, 1:50) with human hTERT-RPE1 (MeOH fixation) cells by Dr. Moshe Kim.



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116441(TSPAN10 Antibody) at dilution of 1:1000