

TSPAN12 antibody



Catalog Number: 116443

Product name

TSPAN12 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSPAN12 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, MCF7 cells. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of TSPAN12 recombinant protein (Accession Number: XM_005250239). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116443(TSPAN12 antibody) at dilution of 1:500