

# UBA2 antibody



Catalog Number: 116467

## Product name

UBA2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

UBA2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Jurkat cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 90 kDa

## Preparation

This antibody was obtained by immunization of UBA2 recombinant protein (Accession Number: NM\_005499). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116467(UBA2 antibody) at dilution of 1:500



IP Result of anti-UBA2 (IP:Catalog No:116467, 4ug; Detection:Catalog No:116467 1:500) with Jurkat cells lysate 3200ug.