

TWIST1 antibody



Catalog Number: 116501

Product name

TWIST1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TWIST1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, MCF7 cells, SGC-7901 cells. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of TWIST1 recombinant protein (Accession Number: XM_011515496). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116501(TWIST1 antibody) at dilution of 1:300