

# UBE2DNL antibody



Catalog Number: 116521

## Product name

UBE2DNL antibody

## Specificity

human; other species not tested.

## Antibody description

UBE2DNL Rabbit Polyclonal antibody. Positive IHC detected in human testis.

## Preparation

This antibody was obtained by immunization of UBE2DNL recombinant protein (Accession Number: BC040290). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human testis using Catalog No:116521(UBE2DNL antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:116521(UBE2DNL antibody) at dilution of 1:50 (under 40x lens)