

# UBE2G1 antibody



Catalog Number: 116525

## Product name

UBE2G1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UBE2G1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of UBE2G1 recombinant protein (Accession Number: BC026288). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116525(UBE2G1 antibody) at dilution of 1:100



Immunofluorescent analysis of HepG2 cells, using UBE2G1 antibody Catalog No:116525 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).