

UBE2K antibody



Catalog Number: 116529

Product name

UBE2K antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBE2K Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse testis tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of UBE2K recombinant protein (Accession Number: NM_005339). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116529(UBE2K antibody) at dilution of 1:400

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116529(UBE2K antibody) at dilution of 1:50 (under 10x lens)