

UBE2O antibody



Catalog Number: 116532

Product name

UBE2O antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

UBE2O Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells, Jurkat cells. Observed molecular weight by Western-blot: 150 kDa, 200-230 kDa

Preparation

This antibody was obtained by immunization of UBE2O recombinant protein (Accession Number: BC007822). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116532(UBE2O antibody) at dilution of 1:500