

# UBE2Q2 antibody

Catalog Number: 116533

## Product name

UBE2Q2 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

## Immunogen

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

UBE2Q2 Rabbit Polyclonal antibody. Positive IHC detected in human skin cancer. Positive WB detected in A431 cells, mouse ovary tissue. Observed molecular weight by Western-blot: 43kd

## Preparation

This antibody was obtained by immunization of UBE2Q2 recombinant protein (Accession Number: NM\_173469). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:116533(UBE2Q2 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human skin cancer using Catalog No:116533(UBE2Q2 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

# UBE2Q2 antibody



Catalog Number: 116533

---

human skin cancer using Catalog

No:116533(UBE2Q2 antibody) at dilution of 1:100  
(under 40x lens)