

# UBE2Z antibody



Catalog Number: 116540

## Product name

UBE2Z antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

UBE2Z Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 45kd

## Preparation

This antibody was obtained by immunization of UBE2Z recombinant protein (Accession Number: BC015169). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116540(UBE2Z antibody) at dilution of 1:300