

USP11 antibody

Catalog Number: 116587

Product name

USP11 antibody

Specificity

Human; other species not tested.

Antibody description

USP11 Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis tissue, human brain tissue, human testis tissue. Positive IF detected in HEK-293 cells. Positive WB detected in A549 cells, HEK-293 cells, human testis tissue, Jurkat cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of USP11 recombinant protein (Accession Number: NM_001371072). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:1000-1:10000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116587(USP11 antibody) at dilution of 1:1000




IP Result of anti-USP11 (IP:Catalog No:116587, 4ug; Detection:Catalog No:116587 1:2000) with HEK-293 cells lysate 2400ug.

USP11 antibody

Catalog Number: 116587



Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using Catalog No:116587(USP11 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HEK-293 cells using Catalog No:116587(USP11 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using Catalog No:116587(USP11 Antibody) at dilution of 1:50 (under 10x lens)

