

# UNC5A antibody

Catalog Number: 116615

## Product name

UNC5A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UNC5A Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF7 cells, PC-3 cells. Positive FC detected in MCF-7 cells. Observed molecular weight by Western-blot: 88-100 kDa

## Preparation

This antibody was obtained by immunization of UNC5A recombinant protein (Accession Number: NM\_133369). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116615(UNC5A antibody) at dilution of 1:1000



1X10<sup>6</sup> MCF-7 cells were stained with 0.2ug UNC5A antibody (Catalog No:116615, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.