

UNC5C-Specific antibody



Catalog Number: 116618

Product name

UNC5C-Specific antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

UNC5C-Specific Rabbit Polyclonal antibody.
Positive WB detected in mouse brain tissue.
Observed molecular weight by Western-blot:
50kd,103kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003728).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116618(UNC5C-Specific antibody) at dilution of 1:1000