

UPK2 antibody



Catalog Number: 116629

Product name

UPK2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UPK2 Rabbit Polyclonal antibody. Positive IHC detected in human bladder tissue. Positive WB detected in mouse testis tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of UPK2 recombinant protein (Accession Number: NM_006760). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116629(UPK2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human bladder using Catalog No:116629(UPK2 antibody) at dilution of 1:50 (under 10x lens)