

# UBL4B antibody



Catalog Number: 116652

## Product name

UBL4B antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

UBL4B Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 22-25kd

## Preparation

This antibody was obtained by immunization of UBL4B recombinant protein (Accession Number: NM\_203412). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116652(UBL4B antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using UBL4B antibody Catalog No:116652 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).