

UCHL5 antibody



Catalog Number: 116673

Product name

UCHL5 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UCHL5 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, mouse brain tissue. Positive WB detected in mouse colon tissue, mouse brain tissue. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of UCHL5 recombinant protein (Accession Number: NM_001199262). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:116673(UCHL5 Antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116673(UCHL5 antibody) at dilution of 1:100 (under 10x lens)